


OCULUS  
SURGICAL

# Super View HTC Lens

Disposable Contact Lens



# Super View HTC Lens


## Hassan-Tornamble Contact Lens, disposable

The familiar clarity and image quality of the Tornambe Contact Lens now with an improved stability mechanism developed by Dr. Tarek Hassan. It has an absorbent customizable ring with four legs to maximize stability, providing the surgeon with superior control and visualization.

### Single-use for crisp, clear view every time

- Complete view of posterior pole – 36° field of view
- Maximum instrument maneuverability due to uniform leg orientation
- No assembly: ready to use
- Excellent stability in surgery through absorbent symmetrical customizable foam ring
- 10 single-use Super View HTC Lenses per box

### Trim to fit



### Excellent image quality




Photo: Prof. Claus Eckardt

To order contact your distributor:

[WWW.OCULUSSURGICAL.COM](http://WWW.OCULUSSURGICAL.COM)

OCULUS Optikgeräte GmbH

Postfach • 35549 Wetzlar • GERMANY

Tel. +49 641 2005-0 • Fax +49 641 2005-299

Email: [support@oculussurgical.de](mailto:support@oculussurgical.de) • [www.oculussurgical.com](http://www.oculussurgical.com)

- OCULUS Surgical USA, [info@oculussurgical.com](mailto:info@oculussurgical.com)
- OCULUS Asia, [info@oculus.hk](mailto:info@oculus.hk)
- OCULUS Brasil, [info@oculusbrasil.com.br](mailto:info@oculusbrasil.com.br)
- OCULUS Czechia, [oculus@oculus.cz](mailto:oculus@oculus.cz)
- OCULUS Iberia, [info@oculus.es](mailto:info@oculus.es)
- OCULUS Turkey, [info@oculus-turkey.com.tr](mailto:info@oculus-turkey.com.tr)



OCULUS is certified by TÜV according to  
DIN EN ISO 13485 MDSAP

